



[research/education]



[construction]



[construction]



[architecture/engineering]



Thomas C. Isabell
Assistant TEM Product Manager
JEOL USA, Inc.

Buildings for Advanced Technology Workshop II
Mesa, Arizona

Analytical Instrumentation Facility Requirements for Nanotechnology

January 21-23, 2004



Outline

Analytical Instrumentation for Nanotechnology

- Examples

Facility Design Requirements

Specifications

Environmental Effects on Instrument Performance

Compensation After Construction

Summary

NANOBUILDINGS.COM

Thomas C. Isabell |




JEOL Core Product Groups

- **Electron Optics**
 - Scanning Electron Microscopes (SEMs)
 - Transmission Electron Microscopes (TEMs)
 - Surface Analysis (Auger, EPMA)
 - Scanned Probe Microscopes (SPMs)
- **Electron Optical Instruments for Semiconductors**
 - Electron Beam Lithography
 - Defect Review Tools
 - Inspection Tools
- **Analytical Instruments**
 - Mass Spectrometers
 - NMRs
 - ESRs



NANOBUILDINGS.COM Thomas C. Isabell | **JEOL**

JSM-7700F Scanning Electron Microscope



- High resolution surface analysis
- Resolution of to 0.6 nm
- Magnification range from 25x to 2,000,000x
- Windows 2000 based operating system

NANOBUILDINGS.COM Thomas C. Isabell | **JEOL**




Combined SEM and Raman



Exterior view of SEM/Raman instrument

SEM image of multi-crystal diamond and local Cathodoluminescence spectrum



Permits both SEM observation and Raman spectroscopy by irradiating the specimen with the electron beam of the SEM and the laser beam of a Raman spectrometer coaxial with the SEM's electron beam

NANOBUILDINGS.COM

Thomas C. Isabell |



JSPM-5200 Scanned Probe Microscope



Can be used in a variety of environments: air, controlled atmosphere, liquids or vacuum

Specimen can be heated or cooled


Variety of STM and AFM modes

NANOBUILDINGS.COM

Thomas C. Isabell |



SPM Example




Direct Observation of DNA
Using an Ultralow
Temperature STM

The double helical structure and
the way the molecules face
each other can clearly be
seen.

Specimen courtesy of
Prof. Tomoji Kawai
Director of Intermaterial Research Center,
Osaka University.

NANOBUILDINGS.COM Thomas C. Isabell | **JEOL**


JBX-9300FS Electron Beam Lithography System



Multiple modes
Automation
Dynamic correction system
eliminates defocusing
resulting from beam
deflection

NANOBUILDINGS.COM Thomas C. Isabell | **JEOL**

InP 40nm Periodic Structure on Calixarene Resist



Specimen courtesy of
Tokyo Institute of Technology
Research Laboratory of Ultra-High Speed
Electronics.


NANOBUILDINGS.COM Thomas C. Isabell | JEOL

JEOL Transmission Electron Microscopes (TEMs)

2100F:	200kV FEG TEM/STEM
2200FS:	200kV FEG with Omega Filter
2500SE:	200kV FEG STEM/TEM
2200FS ACEM:	200kV FEG TEM/STEM with Omega and STEM Cs corrector
2200F Mirai:	200kV FEG TEM/STEM with Omega and Monochromator
2200F AOTA:	200kV FEG TEM/STEM with All of the Above
3200S:	300kV LaB6 TEM with Omega Filter
3200FSC:	300kV FEG TEM with high tilt Liquid He holder

NANOBUILDINGS.COM Thomas C. Isabell | JEOL



JEM-2100F Field Emission TEM



- New computer architecture**
 - Client/server compatible
 - Com or dll interface
 - Windows FLC/Full Capability
- New vibration Isolation**
 - Balanced Center of Gravity
 - New Goniometer Cover
- New Digital STEM**
 - 50MX
 - 0.136nm Resolution on Factory floor during the day
- Piezo Stage built-in**
- Remote Control**
 - all functions

NANOBUILDINGS.COM Thomas C. Isabell | **JEOL**

Discovery of a Nanocone of Boron Nitride


Unlike carbon nanocones, which have a tip angle of 60°, BN cones have tip angle of 39°.

Possibly could be used as a probe for high resolution AFM.

Specimen courtesy of
Dr. Yoshio Bando
 Head of Integrated Strategy Office,
 Integrated Strategy Office,
 National Institute for Materials Science.


NANOBUILDINGS.COM Thomas C. Isabell | **JEOL**

Formation of Xe Nanoparticles




In-situ formation and combination of Xe nanoparticles in aluminum by injection of Xe ions into the Al.

Specimen courtesy of
Dr. Kazuo Furuya
 Supervising Researcher Nanomaterials
 Laboratory,
 National Institute for Materials Science.

NANOBUILDINGS.COM
Thomas C. Isabell | 


Discovery of Fullerene of Boron Nitride





These BN fullerenes were discovered by using a new method in which a high speed electron beam is irradiated at the specimen while being imaged in the TEM.

Where C60 fullerenes have a soccer ball like structure, these BN fullerenes are shaped like a die.

Specimen courtesy of
Dr. Yoshio Bando
 Head of Integrated Strategy Office,
 Integrated Strategy Office,
 National Institute for Materials Science.

NANOBUILDINGS.COM
Thomas C. Isabell | 



Facility Design Requirements

General

- Room Size
- Lighting
- Heating
- Ventilation
- Privacy
- Safety


Special requirements for analytical instrumentation


- Mechanical vibration
- Magnetic fields
- Electrical disturbances
- Environmental changes
 - Temperature
 - Pressure

NANOBUILDINGS.COM Thomas C. Isabell | **JEOL**

Cs Correctors, Monochromators, Omega Filters


JEM-2200FS





NANOBUILDINGS.COM Thomas C. Isabell | **JEOL**

Performance Improvement



After installation After 6 months After 1 year

Improvements from a) to b) are due to improvements to the microscope. Those from b) to c) are due to improvements in the room environment.


David Muller and John Grazul, *Journal of Electron Microscopy* **50(3)** : 219-226 (2001)

NANOBUILDINGS.COM Thomas C. Isabell | JEOL

Environmental Effects Upon Performance

- Vibrations
 - Acoustic noise
 - Mechanical
- Electrical disturbances
- Magnetic fields
 - Electromagnetic interference (EMI)
- Air flow
- Thermal fluctuations
- Pressure fluctuations

NANOBUILDINGS.COM Thomas C. Isabell | JEOL



Acoustic Noise - External


Planes, trains and automobiles
Machinery
Facility noise

- Other laboratory equipment
- People

Room must be acoustically shielded from outside

- Thick walls
- Insulation
- Heavy doors
- Room-in-room design

NANOBUILDINGS.COM Thomas C. Isabell | **JEOL**



Acoustic Noise - Internal

Air conditioning system
Cooling fans
Electronic noise
Human noise

Eliminate as many of the sources of noise as possible

In-room sound dampening system

- Foam
- Baffles

NANOBUILDINGS.COM Thomas C. Isabell | **JEOL**

Improved Stage Stability



New dampers and clamshell

NANOBUILDINGS.COM

Thomas C. Isabell |



Mechanical Vibrations

Large range of resonant frequencies

- Specimen holders tend to resonate at high frequencies
- Microscope lenses resonate at lower (<30Hz) frequencies

Planes, trains and automobiles

Heavy machinery

Elevators

HVAC

Laboratory equipment

Human induced

Associated microscope equipment

- HV supplies
- Cooling supplies
- Vacuum system

NANOBUILDINGS.COM

Thomas C. Isabell |



Mechanical Vibrations

During design, ensure that building frequencies are as high as possible

- Much easier to attenuate high frequency vibrations

Location, location, location

- Away from loading docks, storerooms, stockyards

Isolated pad for facility

- Decouple laboratory from rest of building

Solid foundation for instrument

- Ground floor

Isolated pad for each instrument

- Decouple instrument room from rest of laboratory
- Decouple instrument from auxiliary equipment

Vibration dampening

- Equipment away from instrument
- Instrument itself

NANOBUILDINGS.COM Thomas C. Isabell | JEOL

Electromagnetic Interference

Electricity supplies

- High-voltage transmission lines
- Local substation
 - Transformers

Subways, electrical trains

Three-phase supply

- If loads uneven – current in neutral conductor

Mains

Equipment motors

- Elevators, hoisting cranes, AC plant

Other equipment

- Arc-welders
- Arc-discharge equipment – coaters, etc.

Grounding

Lighting

Monitors

NANOBUILDINGS.COM Thomas C. Isabell | JEOL

Electromagnetic Interference

For a straight conductor, the magnetic field is given by
Ampere's Law:

$$B = 2I/d$$

B = magnetic field in mG
I = current in mA
d = distance in m

Direction of field given by right-hand rule

A 0.3 mG rms magnetic field can be detected in a 0.3 nm STEM image

At 1m distance from a wire carrying 0.5 mA, a 1mG field will be detected – enough to degrade 0.3 nm performance

NANOBUILDINGS.COM
Thomas C. Isabell |

Origin of Ground Currents

(a)

Single phase circuit

(b)


Three phase circuit

Net AC field generated will be proportional to the current δ lost to ground

David Muller and John Grazul, *Journal of Electron Microscopy* **50(3)** : 219-226 (2001)

NANOBUILDINGS.COM
Thomas C. Isabell |

Magnetic Field Generated by Power Bus



Measured magnetic field from a power bus located 3 feet below the floor. The bus was rerouted.


David Muller and John Grazul, *Journal of Electron Microscopy* **50(3)** : 219-226 (2001)

NANOBUILDINGS.COM

Thomas C. Isabell |



Spectrometer Sensitivity



Measured deflection of zero loss peak due to truck at loading dock 20 feet from column.

1 eV deflection seen for each liter of engine capacity.



At 50 feet, spectrometer ceases to function as a truck detector.

David Muller and John Grazul, *Journal of Electron Microscopy* **50(3)** : 219-226 (2001)


NANOBUILDINGS.COM

Thomas C. Isabell |





"Gazebo" Shielding



Gazebo design greatly stabilized beam fluctuations in the 1-10 sec time frame.

Further wrapping of column with bubble wrap or neoprene improves temperature stability.

David Muller and John Grazul, *Journal of Electron Microscopy* **50(3)** : 219-226 (2001)

NANOBUILDINGS.COM

Thomas C. Isabell |



Pressure Fluctuations

Sample under vacuum ($\sim 1 \times 10^{-6}$ Pa) while sample holder is at ambient (~ 105 Pa)

Changes in the room pressure causes a deflection of 0.1 nm/Pa

In typical building, A/C system overpressures by 10 Pa

- Opening and closing of doors changes this


Weather changes cause slow drift

- Typical low pressure cell is 1600 Pa lower than air in front
- This causes 160 nm of drift as the front passes

NANOBUILDINGS.COM

Thomas C. Isabell |





Facility Design – Network


Newest generation of TEMs are all PC driven
 Many use client-server architecture
 Remote microscopy is becoming more important

- Instrument sharing
- Instrument isolation
- Remove operator from microscope environment

Network requirements for remote microscopy

- Fast
 - LAN: 10BaseT, 100BaseT
 - WAN: 1Mbps or more (T1 T2 T3...)
- Secure
- Internet2 is ideal

NANOBUILDINGS.COM Thomas C. Isabell |



Compensations After Construction

- Field cancellation systems**
 - Passive
 - Active
 - Can only correct one point on microscope column
 - Must decide between imaging and spectroscopy
- Vibration dampening systems**
 - Passive
 - Active
 - Both most effective at low frequencies (10Hz and below)
 - Most common vibrations are from AC motors (30Hz, 60Hz and 120Hz)
- Shielding**
 - Mu-metal
- Acoustic dampening**

NANOBUILDINGS.COM Thomas C. Isabell | JEOL



Summary


Facility location

- Close enough to loading dock, but not too close
- Away from airports, train lines, highways
- Away from power lines, power stations

Room location

- Ground floor or basement
- Away from loading dock, storeroom
- Away from elevators
- Away from interfering laboratory equipment
- Decoupled from rest of facility

NANOBUILDINGS.COM Thomas C. Isabell | JEOL



Summary

Room design

- Separate instrument from auxiliary equipment
- Vibration dampening if necessary
- Sound dampening
- No hidden or buried cables
- Minimize airflow over instrument
- Stable temperature
- Stable pressure
- Lighting

NANOBUILDINGS.COM Thomas C. Isabell | JEOL

The (Nearly) Ideal TEM Facility



**Someday
All TEMs Will
Be In Nebraska**

JEOL

NANOBUILDINGS.COM Thomas C. Isabell | **JEOL**